

Sudden Cardiac Arrest Sample Event Viewed In **defibView** Software

| | |
|--|--|
| <p>Event date: 07/13/2004 19:31:07 Event duration: 36 ECG seconds File name: C:\WINDOWS\Temporary Internet Files\OLK3210\DefibView Sample Event_VF_Shock_NSR.dfb</p> | <p>Patient: MRN: SSN:</p> |
|--|--|

EVENT SUMMARY

| | |
|--|--|
| <p>Location: Responder: Original File: DefibView Sample Event_VF_Shock_NSR.dfb Shocks: 1 Patient data: Name: Address: Telephone: Gender: Age: 0 DOB:</p> | |
|--|--|

| | | |
|----------------|----------|------------------------|
| Facility Name: | Printed: | July 13, 2004 18:43:33 |
|----------------|----------|------------------------|


NOT FOR DIAGNOSTIC PURPOSES

Page 1 of 4

| | |
|--|--|
| <p>Event date: 07/13/2004 19:31:07 Event duration: 36 ECG seconds File name: C:\WINDOWS\Temporary Internet Files\OLK3210\DefibView Sample Event_VF_Shock_NSR.dfb</p> | <p>Patient: MRN: SSN:</p> |
|--|--|

NOT FOR DIAGNOSTIC PURPOSES

Scale: mVolts / time




| | | |
|----------------|----------|------------------------|
| Facility Name: | Printed: | July 13, 2004 18:43:33 |
|----------------|----------|------------------------|

NOT FOR DIAGNOSTIC PURPOSES

Page 2 of 4

Sudden Cardiac Arrest Sample Event Viewed In **defibView** Software (cont.)



Event date: 07/13/2004 19:31:07 Patient: _____
 Event duration: 36 ECG seconds MRN: _____
 File name: C:\WINDOWS\Temporary Internet Files\OLK3210\DefibView Sample Event_VF_Shock_NSR.dfb SSN: _____

EVENT TAGS

| Time | Symbol | Text | Data |
|-------------|--------|------------------|------|
| 19:31:07.59 | CONN | Connect Pads | |
| 19:31:07.69 | AN | Analyzing | |
| 19:31:16.82 | SHKADV | Shock Advised | |
| 19:31:16.87 | CH | Charging | |
| 19:31:25.54 | WTS | Waiting to Shock | |
| 19:31:26.08 | SHK | Shock | |
| 19:31:26.08 | SHKDLV | Shock Delivered | |
| 19:31:28.33 | AN | Analyzing | |
| 19:31:39.31 | PAUSE | Paused | |
| 19:31:43.68 | SD | Shutdown | |

Facility Name: _____ Printed: July 13, 2004 18:43:33
 Page 4 of 4

NOT FOR DIAGNOSTIC PURPOSES