


RAP™ | FACE SHIELD

The award-winning sustainability-packaging company RAP have repurposed their design and manufacturing teams to produce a solution that helps protect from airborne transmission. The new RAP™ Face Shield has been designed for a range of uses, from healthcare to personal protection, when in the workplace or out and about. In a world after lockdown, the RAP™ Face Shield will help people adjust confidently to a gradual return to normality.



Manufactured in the E.U.

CASE SIZE: 400

- Conforms with Category 1 PPE Directive EU 2016/425
- Manufactured in the E.U.
- High care, food grade material:
 - BRC AA manufacturing facility
 - Food Standard manufacturing certification
- Lightweight, robust and fully disposable face mask
- Latest anti-mist film - **RAPGUARD 2019** technology
- Easy to fit and wear
- Offers face protection
- Versatile for all business and lifestyle uses

Instructions for RAP™ Face Shield

- Thread the head band through the slots on the shield, **from right to left**.
- Position the shield central to the D slots and then adjust to make the shield curve.
- Adjust the head band to suit your head size.
- This shield should be worn in conjunction with a mouth and nose mask as well as goggles where advised to provide the best protection.

[CLICK HERE TO WATCH STEP-BY-STEP INSTRUCTIONS ON HOW TO ASSEMBLE RAP™ FACE SHIELD](#)

